



Close wound Suction Unit

An effective device for close wound drainage under negative pressure post operatively. Easily compressible large Bellow Chamber with 800 ml capacity, provides long lasting suction.

Provided with option to use single or double catheter.

Available with curved needle of different sizes matching with catheter sizes to meet moderate to heavy drainage needs.

X-Ray Opaque line provided through out the length of the catheter to ascertain its proper placement in the wound.

Packed in double wrapped packing. Sterile ready for use.

The complete unit consists of:

Bellow with connector.

Connecting tube with clamp.

Curve needle with matching catheter.

Spare perforated catheter with 'y' shaped connector.

Item Size	Code
Surgical Drainage Disposable VacuGMS 10Fr	H-VAC10
Surgical Drainage Disposable VacuGMS 12Fr	H-VAC12
Surgical Drainage Disposable VacuGMS 14Fr	H-VAC14
Surgical Drainage Disposable VacuGMS 16Fr	H-VAC16
Surgical Drainage Disposable VacuGMS 18Fr	H-VAC18