

SALIVA EJECTOR



Characteristic	Specification	
Material Composition	Cannula and Filtering Part	PVC (Polyvinyl Chloride), Latex free, Phthalate free
	Wire Insert	Zinc-coated iron wire
Dimensions	Saliva Ejector Length	150 mm
	Cannula External Diameter	6,55 mm
	Tip Length	14,8 mm
	Tip External Diameter	9 mm
Description	Flexible plastic cannula with metal wire insert and fixed filtering part (tip). Cannula has a smooth surface and filtering part has 6 openings.	

Features

- Disposable saliva ejector
- Smooth and flexible tip
- Non-removable tip
- Safe and comfortable
- Easily keep bend
- Provide optimal suction
- Non-toxic
- Non-sterile