

# HEMODIALYSIS BLOOD TUBING SET



#### Medication ports

Large finger guard helps protect the user from accidental needle sticks; venous line features multiple medication ports for optimum versatility and choice in medication delivery.



#### Color-coded clamps

Patient safety is enhanced with clamps that are color-coded for confident identification of the appropriate line.



#### Drip chambers

Drip chambers feature a versatile medication port and, for an extra measure of safety, are anti-foaming to minimize splashing.

## Features

- Sterile, Single use, Non toxic ,pyrogen free and Biocompatible
- Transparent PVC tubing
- Soft PVC blood drip chamber
- Soft PVC pump tubing
- Injection sites for intermittent medication
- Luer lock priming connector for recirculation
- Revolving luerlock connector
- Sterilization Method: Ethylene Oxide