

Blood Transfusion Set



Bulb Latex



Vented Spike



Smooth Flow Controller

Blood Transfusion Set

B.T set is used for transfusion of blood and its components. Blood components must be filtered during transfusion to remove clots and small clumps of platelets and white blood cells that form during collection and storage.

Available with Bulb Latex

Smooth Flow Controller

Spike with Air vent



Features

It have Sharp vented/Non vented spike suitable for Blood bags and all type of standard blood containers. 200 micron blood filter prevents passage of any clot into the chamber.

- Available in 20 Drop
- Transparent & flexible drip chamber (Single or Double)
- Soft & kink resistant PVC tube, 150 cm long
- Choice of Luer lock or Luer slip, with or without needle