

HIGH PRESSURE EXTENSION LINES

Features

- Made from non-toxic non-pyrogenic material
- With or without injection sites
- ETO gas sterilized
- For single use



HIGH PRESSURE EXTENSION LINE



* with Injection Site

Extension Liens Type	Tube Length	Standard Connector	Inner Diameter	Code
High pressure Extension Lines	50 cm	M/F	1 mm or 3 mm or 4 mm or other diameters available upon request	E.L A079000
	100 cm	M/F		E.L A079001
	150 cm	M/F		E.L A079002
	200 cm	M/F		E.L A079003
	300 cm	M/F		E.L A079004

* OTHRE LENGTHS OR CONNECTORS TYPE UPON REQUEST