


# IV Set With Tube Latex



Tube Latex



Chamber



Smooth Flow  
Controller

## IV Set With Tube Latex

It is used to pass the I.V. Fluid directly into the veins. The system of administration employs a drip chamber which prevents air entering into the blood stream and allows estimate and smooth flow of intravenous fluid in our vein.

---

Available with Tube Latex

---

Smooth Flow Controller

---

With and without air vent

---



## Features

It is made up of medical grade non toxic, non- kink able soft P.V.C. The drip set comprises of drip chamber and drip tube. It is sterilized by Ethylene Oxide Gas.

- Luer lock connector
- Individually sterile
- Flow regulator to monitor the flow of I.V. Fluid passing through the drip
- Available in various premium variants