


Nasal Catheter



Nasal Catheter

It consists of twin nasal tips that ensure equal distribution of oxygen in both nostrils. These nasal tips are fixed over soft & kink-resistant loop of PVC tube which is adjustable around the head of the patient. It is connected with a Star Lumen tube.

Connected with 200 cms long Star Lumen tube

Adult, Pediatric

For single patient use



Features

Soft twin prong nasal tips to ensure equal, volume of oxygen to both air passages. It has 99.9% filtration efficiency protection against Bacterial/Viral cross infections.

- It supports low dead Space and low flow resistance.
- Tubing length available: 200 cm
- Available in Sterile and Non Sterile
- Adjustable nasal tips for maximum patient comfort