

# Three Way Stop Cock



6% Luer  
Lock Connection



360° Rotating Tap



Crystal Clear  
Transparent Channel

## Three Way Stop Cock

Three Way Stop Cock facilitates multiple lines through single I.V. access. It has an accurate and continuous flow of channels. It also supports uniform infusion flow without any turbulence or air bubble.

---

Smooth tap rotation at 360 degrees

---

Tested at 3 bar pressure

---

6% Luer Taper for compatibility

---



## Features

It is made up of a transparent medical grade polycarbonate. Can be used with intra-arterial and intra-venous systems.

- Sterile & blister packed
- The product is also available in Lipid Resistant
- Also available with Red, White and Blue color tap & pegs.
- A red tap indicates if used for arterial purposes
- Blue tap can be used for Intra Venous Systems.