

Model : GMS-DC3M



Technical Specifications

- Electric motors for leg, backrest and height adjustable
- Backrest adjustable up to 70°
- Height adjustable from 50 cm - 75 cm
- Hand control unit with intelligent operation
- The chair can be placed in a flat position
- Comfortable mattress with 10 cm thickness covered with leather
- The material of mattress suitable for disinfection by wiping with disinfectant
- Comfortable armrests with movable position
- Fixed footrest
- 4 castors with brake (diameter 125 mm)
- Chair load weight up to 200 kg
- Infusion pole holder
- Dimensions : 190x84 cm approximately
- Mechanical C.P.R

CONTACT US:



Tel.: +203 4877604 / 695
Fax: +203 4863613
P.o. Box : 256 Alex
Tel. & Fax : +202 33022628
Tel. & Fax : +202 33025928



Head Office:
15 Place Saad Zaghloul
Square, 21131 Alexandria, Egypt.
Cairo Branch:
56El Riyadh St., Ground Floor,
Cairo, Egypt.



expgms@yahoo.com.co.uk



www.Ghatwary-medical.com